




Suitical

HAPPY RECOVERY



Protective and covering apparel for pets

**“With our passion for pets,
we love to make their
recovery as comfortable and
stress-free as possible.”**

*Melvin Kok
CEO, Suitical International*

CONTENTS

	Page:
Introduction.....	4
Suitical products.....	6

SUITICAL COVERING APPAREL

- Recovery Suit®.....	8
- Recovery (double) Sleeve(s)®	12
- Comparison Chart	15

SUITICAL COOLING APPAREL

- DRY® Cooling Vest.....	16
--------------------------	----

SUITICAL PRODUCTS:



Original



Patented



Lead/toxins



Animal-friendly



INTRODUCTION

Imagine this: your loyal friend is recovering from skin conditions, hot spots or wounds (first aid, after being spayed, neutered or other surgical procedures). This can be quite stressful already. It would not benefit your pet if you increase that discomfort and stress level even further by restricting their ability to move, eat or drink normally. This is what happens when you use a "cone of shame" (the old fashioned Elizabethan collar).

Hearing clearly is very important for pets. Amplifying the sounds by using a cone can distract their sense of direction and could cause them to bump into something (for example your valuable furniture).

The best way

Here at Suitical, we believe it is better to cover the problem area instead of restricting the head to prevent dogs and cats from getting to their surgical sutures, wounds, hot spots or skin problems.

Child-friendly

Pets are protected from themselves as well as from other animals/children.

Easy to use, Patented Designs

Suitical products acts as a second skin, allowing pets to enjoy their day with little anxiety or disruption. The snug fit combined with the easy to use rear fastening, prevents the pet from wiggling out of the suit or sleeve. The breathable fabric enables air to circulate around the wound for quicker healing, and it fully covers the affected areas to prevent licking, biting or scratching. Covered areas are also kept dry and clean, preventing bacteria and dirt from causing infection.

Only the best, non-toxic materials

When protecting wounds, hot spots or skin conditions it is important to keep it dry and clean. To achieve this, we use the best high quality, organic cotton (94-96 percent) to make the Suit and Sleeve as breathable and lightweight as possible without compromising strength and comfort. For stretch we use the highest quality Invista® Lycra.

**"REDUCE STRESS AND DISCOMFORT;
THIS HELPS FOR A QUICKER RECOVERY"**



SUITICAL™ PRODUCTS



RECOVERY SUIT® DOG AND CAT



The Recovery Suit® for Dogs is a professional alternative to a cone and full body bandage. Suitable for first aid purposes, to protect and cover hot spots, wounds, surgical sutures and bandages, ointments, skin conditions, when in season or incontinence and a lot more...

The design is based on the animal's shape and anatomy with a snug fit result. Manufactured using a breathable stretchy Cotton/Lycra fabric, giving complete freedom of movement without compromising strength or durability. Thanks to the rear closure effortless to fit.

The Recovery Suit® for Cats is different in size and shape compared to the Recovery Suit® for Dogs. For most cats, it can remain fastened when using a litter box, due to the specially tailored space under the tail.

RECOVERY (DOUBLE) SLEEVE(S)®

As in our Recovery Suits, the Recovery Sleeve® is a left/right interchangeable front leg protection for recovering dogs. The Recovery Double Sleeves® protects both front legs. To cover hot spots, wounds, surgical sutures and bandages, ointments, skin conditions and a lot more...

Made from high quality Cotton/Lycra stretch fabric to provide the best comfort possible. Thanks to the adjustable closure around the chest they are effortless to fit. This also applies some pressure around the chest with a snug fit as a result, which helps to reduce the stress level.



DRY® COOLING VEST

Keep dogs Cool by using this comfortable and easy to refill DRY® Cooling vest. After activating the vest by filling it with fresh water, it will immediately start cooling up to 15°C/59°F below ambient air temperature for one to three days. Thanks to the patented system, the vest will never become wet or damp during use. The cooling effect is the result of evaporation and is largely dependent on humidity and airflow.


So, prevent dogs from developing heat stress! During and after walk (in the sun), in the car (not to be used in a parked car), or whenever cooling is required.

Cooling technology by:




PROFESSIONAL COOLING SYSTEM

1. Fill the vest with water;
2. Hot outside air evaporates the water;
3. The evaporation cools waterwick;
4. Surface of the vest cools down;
5. The cooling effect radiates through the membranes of the inner fabric on to your dog's body.



RECOVERY SUIT® DOG AND CAT

Multi-purpose core body protection



94%
Organic Cotton

6%
Invista® Lycra
(stretch)

0%
Lead/toxins

100%
Animal-friendly

"COMFORTABLE RECOVERY WITHOUT ANY COMPROMISES."

The Recovery Suit® is a professional alternative to a cone and full body bandage. Suitable for **first aid** purposes, to cover and protect **hot spots, wounds, surgical sutures** and **bandages, ointments, skin conditions**, when **in season** or **incontinence** and a lot more...


The design is based on the animal's shape and anatomy with a snug fit as result. Manufactured using breathable stretchy Cotton/Lycra fabric, giving complete freedom of movement without compromising strength or durability. Thanks to the rear closure effortless to fit.



Allows use of litter box and a walk outside

When walking a dog, only the rear end of the Suit needs to be opened, simply roll the underside of the Suit inwards and towards the belly and secure with the strap. For most cats, the Suit can remain fastened when using a litter box, due to the specially tailored space under the tail.

- ⊗ Allows pet to move freely in their environment
- ⊗ Pets can eat and drink as normal
- ⊗ Pets are protected from themselves as well as from others
- ⊗ Breathable lightweight four-way stretch fabric
- ⊗ Stress free design enables faster healing
- ⊗ Multi-purpose use
- ⊗ Machine washable
- ⊗ Allows use of litter box
- ⊗ Perfectly snug fit
- ⊗ Exclusive, patented design



FIRST AID, HOT SPOT, WOUND, SUTURE AND SKIN PROTECTION

- Reduces the annoyance and stress your pet may experience, compared to wearing an Elizabethan collar.
- Pets are protected from themselves as well as from others.
- Our exclusive breathable design enables air to circulate around the affected area while still keeping it dry.
- Helps keep patients warm after (and during) surgery.

The inner panel of the Suit is light colored allowing any fluid discharge to be easily seen. If the wound leaks, a porous sheet or incontinence pad can be inserted inside the Recovery Suit® or between the inner panel.

WHEN IN SEASON OR INCONTINENCE

The Recovery Suit® reduces the risk of blood or urine stains in your house. In contrast to standard shorts, it is very difficult for your pet to wriggle out of.

SKIN CONDITIONS

The Recovery Suit® offers protection where it is required. It helps prevent possible damage to the skin caused by scratching or biting.

SIZES AND HOW TO USE

RECOVERY SUIT® - DOG

SIZE	LENGTH (BASE OF NECK TO BASE OF TAIL)	
3XSMALL	22 - 35 cm	8.7 - 13.8 inch
2XSMALL	33 - 42 cm	13.0 - 16.5 inch
XSMALL	40 - 45 cm	15.7 - 17.7 inch
SMALL	43 - 51 cm	16.9 - 20.1 inch
SMALL +	49 - 57 cm	19.3 - 22.4 inch
MEDIUM	55 - 69 cm	21.7 - 27.2 inch
MEDIUM +	61 - 72 cm	24.0 - 28.3 inch
LARGE	67 - 76 cm	26.4 - 29.9 inch
XLARGE	74 - 82 cm	29.1 - 32.3 inch
2XLARGE	80 - 92 cm	31.5 - 36.2 inch




AVAILABLE IN:

-  BLACK
-  BLUE CAMOUFLAGE
-  PINK CAMOUFLAGE

RECOVERY SUIT® - CAT

SIZE	LENGTH (BASE OF NECK TO BASE OF TAIL)	
3XSMALL	29 - 37 cm	11.4 - 14.6 inch
2XSMALL	33 - 42 cm	13.0 - 16.5 inch
XSMALL	40 - 45 cm	15.7 - 17.7 inch
SMALL	43 - 51 cm	16.9 - 20.1 inch



AVAILABLE IN:

-  BLACK CAMOUFLAGE

HOW TO USE:

STEP 1:

With large dogs (long legs) it is easier to insert one front leg followed by the head and then the other leg into the provided openings. With small dogs (short legs), the head must be inserted first into the shirt followed by each front leg (separately) into the openings provided.

STEP 2:

With large dogs (long legs) it is easier to insert one front leg followed by the head and then the other leg into the provided openings. With small dogs (short legs), the head must be inserted first into the shirt followed by each front leg (separately) into the openings provided.

STEP 3:

With large dogs (long legs) it is easier to insert one front leg followed by the head and then the other leg into the provided openings. With small dogs (short legs), the head must be inserted first into the shirt followed by each front leg (separately) into the openings provided.

WHEN NATURE CALLS:

The **Recovery Suit® Dog** is designed for many purposes. It needs to be fully closed when used for light incontinence or after being spayed/neutered. When nature calls, the rear end of the Suit can be opened, simply roll the bottom part of the Suit inwards and towards the belly and secure with the added hold-up strap.

For most cats, the **Recovery Suit® Cat** can remain fastened when using a litter box, due to the specially tailored space under the tail.

"PETS LOVE THEIR FREEDOM,
PLEASE DON'T TAKE THIS AWAY FROM THEM."



RECOVERY (DOUBLE) SLEEVE(S)[®]

Multi-purpose SINGLE or DOUBLE front leg protection



- ⊗ Perfectly snug fit
- ⊗ Exclusive, patented design
- ⊗ Freedom to move as normal
- ⊗ Easy to use
- ⊗ Adjustable

Left / Right exchangeable

97%
Organic Cotton

3%
Invista[®] Lycra
(stretch)

0%
Lead/toxins

100%
Animal-friendly

"REDUCE STRESS, THIS HELPS FOR A QUICKER RECOVERY."

The Recovery Sleeve[®] is an interchangeable front leg protection for recovering dogs. The Recovery Double Sleeves[®] protects both front legs. Suitable to cover and protect hot spots, wounds, surgical sutures and bandages, ointments, skin conditions and a lot more...

Manufactured using breathable stretchy Cotton/Lycra fabric to provide the best comfort. Thanks to the adjustable closure around the chest they are effortless to fit. This also applies some pressure around the chest with a snug fit as a result, which helps to reduce the stress level.


Extra comfort at the elbows

An additional soft protective layer has been added to the Recovery Double Sleeves[®]. This has been designed for dogs that require some additional comfort at the elbows.



SIZES AND HOW TO USE

RECOVERY (DOUBLE) SLEEVE(S)[®]



SIZE	CHEST CIRCUMFERENCE	LENGTH SLEEVE
3XSMALL	27 - 33 cm / 10.6 - 13.0 inch	10 cm / 3.9 inch
2XSMALL	30 - 43 cm / 11.8 - 16.9 inch	12 cm / 4.7 inch
XSMALL	40 - 52 cm / 15.7 - 20.5 inch	17 cm / 6.7 inch
SMALL	45 - 61 cm / 17.7 - 24.0 inch	23 cm / 9.1 inch
MEDIUM	54 - 68 cm / 21.3 - 26.8 inch	26 cm / 10.2 inch
LARGE	64 - 84 cm / 25.2 - 33.1 inch	35 cm / 13.8 inch
XLARGE	80 - 104 cm / 31.5 - 40.9 inch	41 cm - 16.1 inch
2XLARGE	96 - 124 cm / 37.8 - 48.8 inch	48 cm / 18.9 inch

AVAILABLE IN:
● BLACK

HOW TO USE:

STEP 1:

Insert one leg through the sleeve;

STEP 2:

Insert the other leg through the opening on the opposite side;

STEP 3:

Secure around the chest by using the snaps/press buttons (always secure with at least four snaps).

COMPARISON CHART



	Recovery Suit [®] Dog	Recovery Suit [®] Cat	Recovery Double Sleeve(s) [®] Dog	DRY [®] Cooling Vest Dog	Traditional E-Collar
--	--------------------------------	--------------------------------	--	-----------------------------------	----------------------

Used for	Suitical-Products				
Covers the affected area	✓	✓	✓	✓	
Restricts the head					✓
First Aid coverage	✓	✓	✓	✓	
Anxiety management	✓	✓	✓	✓	
Protects and covers: - Hot spot	✓	✓	✓	✓	
- Wound	✓	✓	✓	✓	
- Bandage / IV line	✓	✓	✓	✓	
After surgery:	✓	✓	✓	✓	✓
- When Spayed / Neutered	✓				
- Wound management	✓	✓	✓	✓	✓
Ointment / salve protection	✓	✓	✓		
Allergies / Skin conditions	✓	✓	✓	✓	✓ head only
Incontinence / when in season	✓				
Safe around children	✓	✓	✓	✓	
Safe around pet's best friends	✓	✓	✓	✓	
Puppy/Kitten nursing control	✓				
Comfortable	✓	✓	✓	✓	
Reduces stress	✓	✓	✓	✓	
Amplifies sounds					✓
Freedom to move, eat and drink	✓	✓	✓	✓	
Furniture friendly	✓	✓	✓	✓	

DRY® COOLING VEST

THE ULTIMATE HEAT-STRESS PREVENTION!



Exterior:
50% Polyester Jersey
50% PU Ultra Fresh
Antibacterial



Core:
Inuteq™ technology
Wicking Felt



Interior:
50% Nylon
50% PU



Lead/toxins



Animal-friendly

**“KEEP YOUR DOG COOL AND DRY
THIS HELPS TO AVOID HEAT STRESS.”**

Keep dogs COOL by using this comfortable and easy to refill DRY® Cooling Vest. After activating the vest by filling it with fresh water, it will immediately start cooling up to 15°C/59°F below ambient air temperature for one to three days. Thanks to the patented system, the vest will never become wet or damp during use. The cooling effect is the result of evaporation and is largely dependent on humidity and airflow.

So, prevent dogs from developing heat stress! During and after a walk (in the sun) or whenever cooling is required.



- ❄️ Comfortable heat-stress management!
- ❄️ Cools down up to 15°C/ 59°F below ambient air temperature;
- ❄️ Immediate cooling up to 3 days, 100% DRY° inside and out;
- ❄️ Unique and patented DRY° evaporation cooling technology;
- ❄️ Proven technology by sporters, military and professionals;
- ❄️ No refrigeration or additional gels/chemicals are required;
- ❄️ Antibacterial, machine washable and easy to clean.




Cooling technology by:




PROFESSIONAL COOLING SYSTEM

1. Fill the vest with water;
2. Hot outside air evaporates the water;
3. The evaporation cools waterwick;
4. Surface of the vest cools down;
5. The cooling effect radiates through the membranes of the inner fabric on to your dog's body.



SIZES AND HOW TO USE

DRY® COOLING VEST



SIZE	LENGTH	CHEST CIRCUMFERENCE	NECK CIRCUMFERENCE	WATER REQUIRED
XSMALL	25 - 35 cm / 10 - 14 inch	40 - 60 cm / 16 - 24 inch	< 35 cm / < 14 inch	100 ml / 3.40 oz.
SMALL	35 - 50 cm / 14 - 20 inch	55 - 70 cm / 22 - 28 inch	< 45 cm / < 18 inch	200 ml / 6.75 oz.
MEDIUM	50 - 60 cm / 20 - 24 inch	60 - 80 cm / 24 - 32 inch	< 50 cm / < 20 inch	290 ml / 10.20 oz.
LARGE	60 - 75 cm / 24 - 30 inch	65 - 85 cm / 26 - 34 inch	< 55 cm / < 22 inch	430 ml / 14.50 oz.
XLARGE	70 - 90 cm / 28 - 36 inch	75 - 95 cm / 30 - 38 inch	< 65 cm / < 26 inch	525 ml / 17.70 oz.

AVAILABLE IN:
 **SILVER**

HOW TO ACTIVATE THE DRY® COOLING TECHNOLOGY:



Fill the vest with fresh water
(see size chart)




Distribute the water evenly.



Remove excess water and air then close the valve.



Vest starts to cool immediately for up to three days.



HOW TO USE:

STEP 1:
Insert the head through the opening;

STEP 2:
Secure the vest by closing the buckles from the belly part to the back;

STEP 3:
Enjoy immediate cooling.



Protect, Recover, Care, Comfort

Melvin Kok
CEO

+31 (0) 85 87 69 701
+31 (0) 61 42 16 308
melvin.kok
melvin@suitical.com
www.suitical.com

LeAnn Wagner
United States Representative

+1 (219) 861 2009
LeAnn@suitical.com
www.suitical.com

Visit
Suitical.com



Follow us
on Instagram!

